



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988

LED Surgical Light Elite 7050

Features

- Light head is streamlined, ultra-thin with aerodynamic shape and laminar compatibility.
- Germany made OSRAM LED bulbs and French Gaggione lens guarantee a uniform nature cool lighting in surface and deep cavity procedures.
- Spring arms from ONDAL, Germany which can turn around 360° at all points of articulation without drift, more easier to maneuver and place.
- Revolving outer handle at the light head allow to adjust focus dia and positioning the light head.
- The light have feather touch control panel on the light suspension for adjustment of light intensity, color temperature and on/off.
- Micro computer digital control with shifts luminance 25%-100% for selection

Technical Specification

Technical Specification	Main Head	Satellite Head
Light head dia	70 cm	50 cm
Illumination intensity at 1 m	180,000 Lux	140,000 Lux
Focus dia	15-30 cm	15-30 cm
Focus adjustable	Yes	Yes
Depth of illumination	140 cm	120 cm
Color temperature	4000 K-5000K	4000 K-5000K
CRI (Color Rendering Index)	95	95
Illuminating adjustment	25%-100%	25%-100%
No. of LEDs	60	45
LEDs lifetime	50000 hrs	50000 hrs
Control unit	SMPS	SMPS



CAPITAL MEDICAL EQUIPMENTS SALES & SERVICES

MediworldIndia

contact@mediworldindia.com
mediworld21@gmail.com
+91 98900 67867 or +91 8668334988